

CleanFreak® Endoscope Cleaning Brushes



End to end, bristle for bristle, there's nothing quite like CleanFreak.®

What do you expect from your endoscope cleaning brushes? How about a broad line from one supplier, and confidence in the products that help you clean sensitive, expensive equipment. CleanFreak Endoscope Cleaning Brushes are precisely what you expect from Medovations—solutions at a price that won't clean out your budget. End to end, bristle for bristle, you'll feel better with a CleanFreak by your side.

Heavy Duty Cleaning Brushes

Our three newest brushes now feature a wire reinforced catheter. These heavy-duty brushes have been designed to handle more rigorous use when cleaning endoscopic equipment. Like all CleanFreak® Endoscope Cleaning Brushes, these new heavy-duty solutions feature a protective tip and soft, yet effective bristles that won't harm sensitive surfaces.



Heavy Duty Standard Brush
1200-07
This simple, effective single-head brush combines a non-slip, wire-reinforced catheter with a tapered head and a protective tip. It's the go-to brush for heavy-duty cleaning.



Heavy Duty Dual Ended Brush
1200-08
This heavy-duty dual-ended brush features soft bristles, a protective plastic tip, and a compact head that provides optimal cleaning.



Heavy Duty Combination Brush
1200-09
This convenient dual-ended brush features a valve brush with a comfortable grip on one end and a channel brush on the other. Its wire-reinforced catheter remains firm to glide through longer channels.

Solutions at a price that won't clean out your budget.

You know Medovations for its broad lines of gastroenterology and surgical products. Our line of CleanFreak® brushes is designed to meet the cleaning requirements of flexible endoscopes:

- Range of brush sizes and configurations
- 1.7 mm catheter for control and flexibility
- Single-use, durable design
- Non-sterile and not made with natural rubber latex

CleanFreak® Endoscope Brushes for Flexible Endoscopes

Part No.	Description	Bristle Diameter	Sheath Diameter	Length	Channel Diameter Range	Quantity/Case
1200-01	Standard Brush	5-7 mm	1.7 mm	230 cm	2.0-6.0 mm	50
1200-02	Extra Wide Brush	7-9 mm	1.7 mm	230 cm	2.8-6.0 mm	50
1200-03	Extra Long Brush	5-7 mm	1.7 mm	230 cm	2.0-6.0 mm	50
1200-04	Dual Ended Brush	5-7 mm	1.7 mm	230 cm	2.0-6.0 mm	50
1200-05	Combination Brush (Valve & Channel Brush)	5-7 mm	1.7 mm	230 cm	2.0-6.0 mm	50
1200-06	Valve Brush	1.0 cm				50
1200-07	Heavy Duty Standard Brush	5-7 mm	1.7 mm	230 cm	2.0-6.0 mm	50
1200-08	Heavy Duty Dual Ended Brush	5-7 mm	1.7 mm	230 cm	2.0-6.0 mm	50
1200-09	Heavy Duty Combination Brush (Valve & Channel Brush)	5-7 mm	1.7 mm	230 cm	2.0-6.0 mm	50



Standard Brush 1200-01

This single-head brush is a must have accessory, combining a non-slip catheter with a tapered head and a protective tip. If one brush ever did the work of many, this would be it.



Extra Wide Brush 1200-02

This brush contains longer bristles, a plastic tip, and a compact head that provides optimal cleaning capability for channel diameters between 2.8 mm and 6 mm.



Extra Long Brush 1200-03

The tapered brush head and firm yellow catheter flexes easily when guided through the channel. With a longer head design, this brush doubles the cleaning contact with every pass through the channel.



Dual Ended Brush 1200-04

This dual-ended brush enables efficient cleaning by using two brushes in one pass. The textured catheter helps maintain a secure grip.



Combination Brush 1200-05

For maximum convenience, a standard brush and a valve brush are now in one product. A bright yellow tip on the channel brush helps protect the channel from abrasion.




Valve Brush 1200-06

This compact cleaning brush features a densely packed head of soft, gentle bristles to clean the delicate areas around the control head and valve cavities. The handle design ensures a comfortable, secure grip.



Esophageal Dilators



CleanFreak® Cleaning Products



Endoscope Care & Accessories

At Medovations, we are committed to providing you with high-quality, durable products for your gastrointestinal procedures. See the complete line of GI solutions at **Medovations.com**

Learn more from your personal product representative at **800-558-6408** or **Medovations.com**



Medovations, Inc. | 102 East Keefe Avenue | Milwaukee, Wisconsin 53212 USA
800.558.6408 or 414.265.7620 | Fax 414.265.7628
Orders: **orders@medovations.com** | General Inquiries: **medo@medovations.com**