



Universal Funnel Tubing



- Clear, kink-resistant PVC
- "Bubble" every 3ft (1m) provides QuickConnect™ integral connector for devices, suction units, or tubing
- Packaged in 100ft (30m) roll in an easy-dispensing storage box
- High resistance to vacuum collapse



"Bubble" every 3ft (1m)



Integral connector formed when "Bubble" is cut

Part Number	I.D.	Length	Qty per cs
1268-01	3/16"(5mm)	100'(30m)	1 Roll
1268-02	1/4"(6mm)	100'(30m)	1 Roll
1268-03	9/32"(7mm)	100'(30m)	1 Roll